


## violetFluor™ Antibodies

Find success with violetFluor™ reagents! These violet dye antibodies offer an economical solution for your violet laser channels. With performance and brightness equivalent to top alternatives, violetFluor™ antibodies give you the brightness you need without breaking your budget.




*Mouse spleen lymphocytes were stained with Anti-Mouse CD4 (RM4-5) using violetFluor™ 500 (cat. no. 85-0042) (left) and a top alternative (right).*

## violetFluor™ 450 Reagents

Excitation: 405 nm | Emission: 450 nm

Description	Clone	Sizes	Cat. No.
Anti-Human CD2	RPA-2.10	25 tests, 100 tests, 500 tests	75-0029
Anti-Human CD3	OKT3	25 tests, 100 tests	75-0037
Anti-Human CD3	SK7	25 tests, 100 tests	75-0036
Anti-Human CD3	UCHT1	25 tests, 100 tests	75-0038
Anti-Human CD4	OKT4	25 tests, 100 tests	75-0048
Anti-Human CD4	RPA-T4	25 tests, 100 tests	75-0049
Anti-Human CD4	SK3	25 tests, 100 tests	75-0047
Anti-Human CD8	SK1	25 tests, 100 tests, 500 tests	75-0087
Anti-Human CD8a	Hit8a	25 tests, 100 tests	75-0089
Anti-Human CD8a	OKT8	25 tests, 100 tests	75-0086
Anti-Human CD8a	RPA-T8	25 tests, 100 tests	75-0088
Anti-Human CD11c	3.9	25 tests, 100 tests	75-0116
Anti-Human CD14	61D3	25 tests, 100 tests	75-0149
Anti-Human CD16	3G8	25 tests, 100 tests	75-0166
Anti-Human CD19	HIB19	25 tests, 100 tests	75-0199
Anti-Human CD20	2H7	25 tests, 100 tests, 500 tests	75-0209



## violetFluor™ 450 Reagents (continued)

Description	Clone	Sizes	Cat. No.
Anti-Human CD27	O323	25 tests, 100 tests	75-0279
Anti-Human CD38	HIT2	25 tests, 100 tests	75-0389
Anti-Human CD45	HI30	25 tests, 100 tests	75-0459
Anti-Human CD45RA	HI100	25 tests, 100 tests	75-0458
Anti-Human CD56 (NCAM)	MY31	25 tests, 100 tests	75-0564
Anti-Human CD161	HP-3G10	25 tests, 100 tests	75-1619

View complete list at: [tonbobio.com/collections/violetfluor-antibodies](https://tonbobio.com/collections/violetfluor-antibodies)

## violetFluor™ 500 Reagents

Excitation: 405 nm | Emission: 500 nm

Description	Clone	Sizes	Cat. No.
Anti-Human CD3	OKT3	25 tests, 100 tests	85-0037
Anti-Human CD3	UCHT1	25 tests, 100 tests	85-0038
Anti-Human CD4	OKT4	25 tests, 100 tests	85-0048
Anti-Human CD4	RPA-T4	25 tests, 100 tests	85-0049
Anti-Human CD4	SK3	25 tests, 100 tests	85-0047
Anti-Human CD8	RPA-T8	25 tests, 100 tests	85-0088
Anti-Human CD8a	OKT8	25 tests, 100 tests	85-0086
Anti-Human CD14	61D3	25 tests, 100 tests	85-0149
Anti-Human CD16	3G8	25 tests, 100 tests	85-0166
Anti-Human CD20	2H7	25 tests, 100 tests	85-0209
Anti-Human CD45	HI30	25 tests, 100 tests	85-0459
Anti-Human/Mouse CD45R (B220)	RA3-6B2	25 µg, 100 µg	85-0452
Anti-Mouse CD3	17A2	25 µg, 100 µg	85-0032
Anti-Mouse CD4	GK1.5	25 µg, 100 µg	85-0041
Anti-Mouse CD4	RM4-5	25 µg, 100 µg	85-0042
Anti-Mouse CD8a	53-6.7	25 µg, 100 µg	85-0081
Anti-Mouse CD8a	2.43	25 µg, 100 µg	85-1886
Anti-Mouse CD11b	M1/70	25 µg, 100 µg	85-0112
Anti-Mouse CD44	IM7	25 µg, 100 µg	85-0441
Anti-Mouse CD45	30-F11	25 µg, 100 µg	85-0451
Anti-Mouse CD45.1	A20	25 µg, 100 µg	85-0453
Anti-Mouse CD45.2	104	25 µg, 100 µg	85-0454
Anti-Mouse Ly-6G	1A8	25 µg, 100 µg	85-1276

View complete list at: [tonbobio.com/collections/violetfluor-antibodies](https://tonbobio.com/collections/violetfluor-antibodies)

